





CEILING LIFT SYSTEM


MODEL: MD-T-200

SHANGHAI MEDDO MEDICAL DEVICE CO., LTD

 meddomedical2012@gmail.com

 We chat:+8617715168788

 +86-17715168788

 www.meddomedical.com



Shanghai Meddo Medical Devices Co., Ltd. (Meddo), established in 2012, is a high-tech enterprise dedicated to R&D, sales and technical service in clinical nursing, emergency medical aid, infection control in wards and medical rehabilitation such as Patient lift system, Automatic Cleaning and Disinfection Machine-Bedpan Washer etc.

We provide the customers with OEM&ODM customized service based on the cost-effective solutions, and will make full use of our R&D technical team and cost control force to provide overseas customers with high-quality and cost-effective medical and maintenance products, infection control equipment for wards, critical first-aid and rehabilitation equipment and intelligent medical equipment, benefiting the world with high-quality medical products and services.

Smart Ceiling Lift system MD-T-200

MD-T-200 is specifically engineered for homes, hospitals, and professional care facilities, setting a new standard for residential and institutional care. Designed for daily sitting and lying transfers, it delivers practical support for both in-home care and clinical care scenarios.

MD-T-200 is centered around a lightweight and multifunctional overhead track, utilizing a modular design that allows for quick installation, eliminating the need for complex engineering modifications. The track system can be flexibly arranged according to space requirements, ensuring efficient coverage without dead zones within the room. Whether it's from the hospital bed to the bathroom, from the wheelchair to the sofa, or from indoors to the hallway, it provides smooth and continuous transfer support, greatly optimizing the nursing workflow and enhancing care efficiency.



■ Comfortable Experience

Smooth angle adjustment, providing full care for the user's posture and dignity throughout the process.

■ Effective Nursing

Quick installation and intuitive operation significantly reduce the burden of caregiving.

■ Safe and Reliable

Electric power assistance paired with professional slings provides dual safeguards for safe transfer.

■ Flexible Switching

Fixed/portable dual mode one device for multiple scenarios.

■ All-domain Suitable

Hospitals, nursing homes, and home care services provide comprehensive coverage.

Safe Transfers, Total Peace of Mind

Extremely lightweight and compact. Focuses on sitting and lying transfer care

Whether moving from the bedroom to the bathroom, living room to balcony at home, or across different zones in hospitals and care institutions. MD-T-200 ensures smooth and secure transfers for sitting and lying positions across all indoor areas. It reduces the risk of secondary injuries such as falls or bumps caused by manual lifting, ensuring safety throughout the process.



Transfer in Home



Transfer in Hospital




Toileting Support



Training Support

Unique Innovation

One machine, dual-purpose, switching between fixed and portable at will



MD-T-200 delivers remarkable functional innovation: the same equipment unit can serve as a fixed overhead track system, and also be easily converted into a portable mobile system.

By simply flipping the main unit, users can freely switch between upright hanging and reverse hanging modes. This design breaks scene constraints and frees care from spatial limitations. Whether moving between multiple rooms at home or between different care settings, it provides continuous, stable and professional support.


Humanized Design

Multi-angle adjustment, ensuring comfort throughout the journey

Our system features a lightweight, multifunctional overhead track with modular design for quick installation, no complex modifications needed. The track can be flexibly arranged to ensure full room coverage without dead zones. It enables smooth, continuous transfers between beds, bathrooms, wheelchairs, sofas and hallways, streamlining nursing workflows and improving care efficiency.



Example for Home Design



Barrier-free Restroom

Enhanced Quality of Life, Preserved Dignity

Greatly lowering the subsequent use burden for both families and institutions.



The design of the portable support rod fully considers ergonomics and psychological needs: it can automatically maintain the optimal support angle throughout the transfer process, ensuring comfortable sitting posture and even support for the user's body, effectively alleviating local pressure. At the same time, this design creates a natural eye level for caregivers, allowing them to easily maintain eye contact with the user during the operation, conveying comfort and encouragement, thereby enhancing mutual trust and making the care process more warm and respectful.

A Variety of Slings

Safe and Comfortable ,Support Diversified Scenarios


The two-point hanger can be used with slings (including walking training patient slings and sitting patient slings, Supine Sling) for showering, toileting, and prone position ventilation.



For example, assisting users with toileting, showering, wheelchair transfer, flat transfer, and prone position ventilation treatment improves the work efficiency of caregivers, protects their lumbar vertebrae, and avoids possible secondary injuries to patients during movement.

Various Fittings Support

A Wide Range of Applications



The door passing system and layout scheme can be customized to make it compatible with different facilities (fire protection systems, air conditioners, lights, and doors), meeting the needs of various occasions.




Smart Ceiling Lift system

MD-T-200



Specification:

- Maximum Lifting Weight: 200kg
- Belt Length: Max 2150 mm
- Weight: 3.8kg (without pole)
- Size: 198 x 188 x 98mm
- Waterproof Rating: IPx4
- Warranty: 5 years



For more information, please contact our professional technical staff.