

Sponge Mattress

Product Details-Sponge Mattress

Shanghai Meddo Medical Devices Co.,Ltd





Features








- TDI BASF and polyurethane foam sponge, prevents patients from developing pressure sore
- Maintenance-free, 5 years lifetime
- Fire resistant comply with CA117/BS7177
- Maintains air flow to reduce heat and moisture.

- Water Retarcant
- Environmental friendly
- Customizable







Technical specification

Cover

-  · PU Material, Two-stretch elastic fabric
-  · Waterproof surface coating, Vapour Permeable
-  · Antibacterial, Antimicrobial, bacteriostatic and fungistatic tested
-  · Hidden Zipper
-  · Washable

Foam mattress

-  · With the x-types fiber designs, double-sides mesh; present a 6-sides breathable 3D structure
-  · Breathable
-  · Effectively releases cervical pressure, and also provide a great support to the body. Keeps the skin in dry
-  · Washable

Shanghai Meddo Medical Devices Co.,Ltd

Tel:17715168788 Website:www.meddomedical.com Linkenin:MeddoMedical