

Circuits for Medijet[®] and Nasal Canulas *



REF	Name	Diameter Connectors	Inspiration/expiration	Additional Connectors	Heated	Accessories	Single use/reuse
1207	Patient circuit	F22,M7.4	Inspiration	1x softconnector; 1x silicone piece M22/ M15; F22/ M22;M22/M22	Yse-for F&P 850	No	Single use
1207 MKI	Patient circuit	F22,M7.4	Inspiration	1x softconnector; 1x silicone piece M22/ M15; F22/ M22;M22/M22	Yse-for F&P 850	Chamber	Single use



Medijet[®] *

REF	Name	Single use/reuse
1000	Medijet [®] (20pieces/unit)	Single use
1000-10	Medijet [®] Starter Set(for 10 patients)	Single use
1000-20	Medijet [®] Starter Set(for 20 patients)	Single use



Prongs *

Size	x-small	small	medium	large	x-large	medium wide	large wide
REF	1200-01	1200-021	1200-02	1200-022	1200-03	1200-032	1200-033
Nostril Ø [mm]	3	3.5	4.1	4.75	5.5	3.7	5
Packing unit	10	10	10	10	10	10	10



Masks *

Size	small	medium	large	x-large
REF	1200-04	1200-05	1200-06	1200-07
Mask width nose [mm]	18.5	21.0	23.0	26.0
Packing unit	10	10	10	5



Bonnets-Disposable *

*Made by Medin

Size	XXS	XS	S	M	L	XL	XXL	XXXL
REF	1213-10	1214-10	1215-10	1216-10	1217-10	1218-10	1219-10	1220-10
Head circumference [mm]	17-19	19-21	21-23	23-25.5	25.5-28	28-30	30-33	33-36
Packing unit	10	10	10	10	10	10	10	10



MEDDO-CPAP

Comprehensive respiratory support for pretermatures and neonates from the firstminute of life



For more information, please contact us

177-1516-8788

SHANGHAI MEDDO MEDICAL DEVICES CO.,LTD







Linkenin:MeddoMedical

www.meddomedical.com




Overall Excellent Breathing Technique


The ventilator is the basic NCPAP device, It has NCPAP(with Medijet) ventilation mode and provides a reliable, economical solution for basic NIV therapy.

-  Spontaneous respiratory rate display
-  Low noise
-  5 hours of operation on battery power
-  Pressure and oxygen measurement with alarm
-  high-flow oxygen therapy
-  Integrated O2/air blender

MEDDO-CPAP Mode & Specification



NCPAP Mode



High-flow Oxygen Therapy

Specification	
Size	18*24*14.5cm
Weight	2.6kg
Gas Delivery Unit	Mechanical Blender
Gas Supply	3.5-6.0 bar (Air and oxygen)
Gas Connection	DISS/NIST
Power	Internal Power Supply, 9V DC, 5 hours Working time, Rechargeable
	External Power Supply, 100-240V AC 50-60Hz/ 400mA
	Mechanical Blender provides backup support, constant flow, no alarm
Display	Monochrome, 96*61mm
CPAP	0-10mbar
Alarm	1-99bpm; In CPAP mode
Safety	Mechanical overpressure valve(80cmH2O)
Date Interface	RS232
Tendency	FiO2